



Portable Coffee Refractometer LD-LCR-A10

Labodam Equipment Ltd.

info@labodam.com | www.labodam.com

Portable Coffee Refractometer LD-LCR-A10

Overview

Portable coffee refractometer LCR-A10 is a microprocessor controlled pocket sized optical instrument that provides instant measurement readings. LCD screen displays measurement data in one view. Automatically turns off in 1 minute of non-use prolonging it's battery life. It is designed for measuring brix and TDS of coffee.

Features:

- Digital meter displays easy to read results
- Standardized Brix scale for easy read and accuracy
- Automatic and zero calibration with distilled or tap water
- Abrasion and corrosion-resistant coating on sapphire prism
- Auto power off (after 1 minute)

Applications:

Used in laboratories, research institutes, beverage industry, food industry, etc.

Specifications:

Weight	90 grams
Battery	1.5 V
HS Code	90275000
TDS range	0.00 ~ 25.00
Brix range	0.00 ~ 30 %
Dimensions	121 × 58 × 25 mm
Temperature	10 ? ~ 100 ?