



Vertical Laboratory Autoclave LD-LVA-K13

Labodam Equipment Ltd.

info@labodam.com | www.labodam.com

Vertical Laboratory Autoclave LD-LVA-K13

Overview

Vertical Autoclave LVA-K13 is a top loading sterilizer which provides safe, economic and effective sterilization of laboratory instruments. The autoclave has automatic system for feeding water, pulse vacuum, sterilization, vacuum drying, exhaustion and draining. These autoclave comes with digital display with touch key pad.

Features:

- 7 inches' touch screen display
- Automatic control system
- Elevated type door
- Automatic system for feeding water, pulse vacuum, sterilization, vacuum drying, exhaustion and draining
- Equipped with PT-TT testing port
- Energy efficient steam generator
- Equipped with 0.2µm sterile ventilation filter
- Clean and dry sterilizing environment facilitated by steam-water inner circulation system

Applications:

This product finds applications in: Biochemistry, Pharmaceutical Industry, Medicine and Health care Industry, Biotechnology, Microbiology and Scientific Research.

Specifications:

Capacity	120 L
Timer Range	0 - 99 min
Heat Average	± 1
Sterilization Pressure	0.22 MPa
Sterilization Temperature	134 °C
Chamber Dimension (H × D)	504 x 404 x 590 mm