



Floor Balance LD-LFRB-A20

Labodam Equipment Ltd.

info@labodam.com | www.labodam.com

Floor Balance LD-LFRB-A20

Overview

Floor Balance LD-LFRB-A20 is an electronic balance comprising indicator and square weighing pan, offers 50 kg maximum weighing capacity. Features 304 stainless steel iron pan with ABS house and aluminum sensor, detachable indicator provides kg, g, oz, lb as weighing units. Equipped with LCD screen display and RS232 interface, this equipment offers counting, percentage, full range tare, animal weighing function with over-load protection, zero tracking and low-battery indicator. With quick response, accurate weighing and high quality, this equipment brings great convenient to user.

Features:

- Electronic balance with detachable indicator and square weighing pan
- 304 stainless steel iron pan with ABS house and aluminum sensor
- Weighing Units: kg, g, oz, lb
- LCD screen display with white backlight and RS232 interface
- Offers counting, percentage, full range tare, animal weighing function
- Over-load protection, zero tracking and low-battery indicator
- Quick response, accurate weighing and high quality
- Memory of last setting for convenience

Applications:

Floor electronic balances are used in the food, pharmaceutical, textile, geological and jewellery industries, as well as in schools, hospitals, laboratories, medical, chemical and other large factories.

Specifications:

Power	AC
Speed	1 s
Sensor	Aluminum sensor
Display	LCD screen display
Accuracy	1 g
Capacity	50 kg
Material	304 Stainless steel Iron, ABS house
Pan type	Square
Functions	Counting, percentage, full range tare, animal weighing function, zero tracking
Indicator	Low battery
Net weight	6 kg

Gross weight	7 kg
Weighing units	kg, g, oz, lb
Value stable time	2 s
Overload protection	Yes
Pan to floor height	120 mm
Pan dimensions (L×B)	400×500 mm
Indicator dimensions (L×H)	250×150 mm
Packaging dimensions (W×D×H)	490×450×170 mm