



Floor Balance LD-LFRB-A21

Labodam Equipment Ltd.
info@labodam.com | www.labodam.com

Floor Balance LD-LFRB-A21

Overview

Floor Balance LD-LFRB-A21 is made to offer reliable and practical weighing for industrial applications, quality control settings, and labs. Routine weighing is made easy and effective by its useful design, which facilitates efficient daily laboratory and manufacturing workflows.

Features:

- Accurate and reliable weighing performance
- Fast response for efficient routine weighing
- Clear and easy-to-read backlit LCD display
- Multiple weighing-unit selection for greater flexibility
- Unit-status memory function for convenient operation
- Calibration-value selection
- Automatic zero-tracking for stable measurements
- Counting function for quantity-based weighing
- Percentage weighing function

Applications:

Suitable for routine weighing of samples, materials, and containers in laboratories, quality control, manufacturing, and industrial environments.

Specifications:

Pan Size	30 cm × 40 cm
Interface	RS232 – Optional Accessory
Readability	0.1 g
Response Time	? 3 s
Linearity Error	? 0.3 g
Maximum Capacity	30 kg
Corner Load Error	? 0.3 g