



## **Rapid Low Temperature Circulating Bath LD-LRTB-A11**

**Labodam Equipment Ltd.**

**[info@labodam.com](mailto:info@labodam.com) | [www.labodam.com](http://www.labodam.com)**

# Rapid Low Temperature Circulating Bath LD-LRTB-A11

## Overview

Rapid Low Temperature Circulating Bath LRTB-A11 is a microcomputer controlled 3 L bath with cold ring liquid circulation and cryogenic liquid with a temperature range of  $-10\text{ }^{\circ}\text{C}$ . Refrigeration unit with low temperature cooling pump facilitates consistent high pressure and maintenance of low temperature under working conditions of reaction experiments. Equipped with anti-corrosive circulatory system preventing traffic flow, high pressure and low pump body temperature.

## Features:

- Volume - 3 L
- Low temperature cooling pump
- Microcomputer controlled temperature
- Temperature and time LCD display
- Built-in liquid level monitor
- Refrigeration unit with relay , protector and capacitor
- Anti-corrosive circulatory system

## Specifications:

|                        |                                   |
|------------------------|-----------------------------------|
| Volume                 | 3 L                               |
| Pump flow              | 20 L / min                        |
| Head of delivery       | 10 m                              |
| Temperature range      | - 10 °C                           |
| Temperature stability  | $\pm 0.5\text{ }^{\circ}\text{C}$ |
| Refrigerating capacity | 673 W                             |