



Rapid Low Temperature Circulating Bath LD-LRTB-A40

Labodam Equipment Ltd.

info@labodam.com | www.labodam.com

Rapid Low Temperature Circulating Bath LD-LRTB-A40

Overview

Rapid Low Temperature Circulating Bath LRTB-A40 is a microcomputer controlled 7 L bath with refrigeration system and a temperature range of - 40 °C to 200 °C. Equipped with dual red and green LED windows with upper window measuring red value and lower measuring green value. Features internal and external cycles with external cycles discharging thermostatic liquid outside the bath raising the temperature of the external vessel.

Features:

- Volume - 7 L
- Low temperature cooling pump
- Microcomputer controlled temperature
- Temperature and time LCD display
- Built-in liquid level monitor
- Refrigeration unit with relay , protector and capacitor
- Anti-corrosive circulatory system

Specifications:

Volume	7 L
Pump flow	20 L / min
Heating power	1200 W
Head of delivery	10 m
Temperature range	- 40 °C to 200 °C
Temperature stability	± 0.1 °C
Refrigerating capacity	2329 W
Temperature control tier number	32