



Rapid Low Temperature Circulating Bath LD-LRTB-A70

Labodam Equipment Ltd.

info@labodam.com | www.labodam.com

Rapid Low Temperature Circulating Bath LD-LRTB-A70

Overview

Rapid Low Temperature Circulating Bath LRTB-A70 is a microcomputer controlled 20 L bath with cold ring liquid circulation and cryogenic liquid with a temperature range of - 10 °C. Refrigeration unit with low temperature cooling pump facilitates consistent high pressure and maintenance of low temperature under working conditions of reaction experiments. Equipped with anti-corrosive circulatory system preventing traffic flow, high pressure and low pump body temperature.

Features:

- Volume - 20 L
- Low temperature cooling pump
- Microcomputer controlled temperature
- Temperature and time LCD display
- Built-in liquid level monitor
- Refrigeration unit with relay , protector and capacitor
- Anti-corrosive circulatory system

Specifications:

Volume	20 L
Pump flow	20 L / min
Head of delivery	10 m
Temperature range	- 10 °C
Temperature stability	± 0.5 °C
Refrigerating capacity	3960 W