



## **Rapid Low Temperature Circulating Bath LD-LRTB-A72**

**Labodam Equipment Ltd.**

**[info@labodam.com](mailto:info@labodam.com) | [www.labodam.com](http://www.labodam.com)**

# Rapid Low Temperature Circulating Bath LD-LRTB-A72

## Overview

Rapid Low Temperature Circulating Bath LRTB-A72 is a microcomputer controlled 20 L bath with cold ring liquid circulation and cryogenic liquid with a temperature range of - 40 °C. Refrigeration unit with low temperature cooling pump facilitates consistent high pressure and maintenance of low temperature under working conditions of reaction experiments. Equipped with anti-corrosive circulatory system preventing traffic flow, high pressure and low pump body temperature.

## Features:

- Volume - 20 L
- Low temperature cooling pump
- Microcomputer controlled temperature
- Temperature and time LCD display
- Built-in liquid level monitor
- Refrigeration unit with relay , protector and capacitor
- Anti-corrosive circulatory system

## Specifications:

|                        |            |
|------------------------|------------|
| Volume                 | 20 L       |
| Pump flow              | 20 L / min |
| Head of delivery       | 10 m       |
| Temperature range      | - 40 °C    |
| Temperature stability  | ± 0.5 °C   |
| Refrigerating capacity | 7939 W     |