



High speed centrifuge LD-LHS-B21

Labodam Equipment Ltd.

info@labodam.com | www.labodam.com

High speed centrifuge LD-LHS-B21

Overview

High speed centrifuge LD-LHS-B21 is a versatile device for wide range of application, this centrifuge is provided to give more throughput and greater control. It can accommodate various types rotors according to user. The brushless DC motor enables smooth operation with low maintenance. High speed centrifuge is useful in isolation of sample that require higher range of centrifugal force for their separation.

Features:

- LED display for observation of parameters such as speed, time and RCF
- Compact design low noise operation
- Various kinds of rotor available
- Automatic electric lid lock and speed imbalance protection
- High grade of steel for safe and reliable operation

Applications:

It is widely used in protein precipitation, RNA extraction, Clinical sample processing.

Specifications:

Noise	≤60 dBA
Timer	0~99 min
Max RCF	30910 x g
Dimension	700 x 500 x 400 mm
Max Speed	21000 rpm
Max Volume	4 x 600 ml
Net Weight	67 kgs
Power supply	AC 220 / 110 V 50 / 60Hz
Speed Accuracy	±20 rpm