



Solar freezer LD-LSF-A11

Labodam Equipment Ltd.

info@labodam.com | www.labodam.com

Solar freezer LD-LSF-A11

Overview

This device is fast cooling and high temperature resistant, with advance design, low maintenance for quality performance, safety and accurate results.

Features:

- Fast cooling with optional temperature control
- Freezer is equipped with drawers for safe storage
- Automatic low voltage disconnection for protection
- Low maintenance and easy to clean
- Design with single top door with door lock to prevent unauthorized access
- Environment friendly CFC free refrigerant R134a is used
- Designed to ensure high temperature resistance and minimum energy consumption
- Equipped with hermetic domus compressor with integrated electronic control panel

Applications:

Widely used in body fluid banks, food industries, life science, biotechnology, genetics, research and development.

Specifications:

Capacity	150 Litres
Net Weight	35 kg
Carton Size	580 × 610 × 870 mm
Gross Weight	38 kg
Input Voltage	DC 12 V / 24 V
Additional Parts	Solar panel: 200 W; Battery: 150 Ah; Charge controller: 12 V / 24 V - 15 A
Power Consumption	75 W
Product Size (L × W × H)	738 × 535 × 845 mm