



Solar freezer LD-LSF-A14

Labodam Equipment Ltd.

info@labodam.com | www.labodam.com

Solar freezer LD-LSF-A14

Overview

This device is fast cooling and high temperature resistant, with advance design, low maintenance for quality performance, safety and accurate results.

Features:

- Environment friendly CFC free refrigerant R134a is used
- Design with single top door with door lock to prevent unauthorized access
- Automatic low voltage disconnection for protection
- Low maintenance and easy to clean
- Fast cooling with optional temperature control
- Freezer is equipped with drawers for safe storage
- Designed to ensure high temperature resistance and minimum energy consumption
- Equipped with hermetic compressor with integrated electronic control panel

Applications:

Widely used in body fluid banks, food industries, life science, biotechnology, genetics, research and development.

Specifications:

Capacity	358 Litres
Net Weight	48 kg
Carton Size	1344 × 610 × 870 mm
Gross Weight	51 kg
Input Voltage	DC 12 V / 24 V
Additional Parts	Solar panel: 200 W; Battery: 150 Ah; Charge controller: 12 V / 24V - 15 A
Power Consumption	120 W
Product Size (L × W × H)	1290 × 598 × 850 mm