



Mini Handheld Homogenizer LD-LMHH-A10

Labodam Equipment Ltd.

info@labodam.com | www.labodam.com

Mini Handheld Homogenizer LD-LMHH-A10

Overview

Mini Hand-held Homogenizer LX140MHH is a hand-held homogenizer with an easy to dismantle grinding pestle. Equipped with adjustable speed for user control.

Features:

- Adjustable speed
- Efficient and effective micro-DC motor for grinding trace tissues
- Reuseable grinding pestle with convenient dismantling

Applications:

Used for grinding various plant tissue sample such as root, stem, leaf, flower, fruit, seed, etc and for grinding various animal tissue sample such as cerebrum, heart, lung, stomach, liver, thymus, kidney, intestine, lymph nodules, muscle, skeleton, etc. Suitable for homogenate.

Specifications:

Weight	0.18 kg
Battery	AA battery x 2
Dimensions	143 x 25 x 143 mm
Speed range	18000 rpm
Working noise	55 dB
Battery charging time	> 2 hours
Battery non-loaded working time	2 hours