



Multi-tube Vortex Mixer LD-LMVM-A12

Labodam Equipment Ltd.
info@labodam.com | www.labodam.com

Multi-tube Vortex Mixer LD-LMVM-A12

Overview

Multi-tube Vortex Mixer LD-LMVM-A12 has a speed range of 500 rpm to 3000 rpm. It uses a brushless DC motor for reliable and smooth operation. The mixer offers touch mode and continuous mode for operation and features an LED display panel for precise control of speed. Its intelligent microprocessor control ensures stable and accurate operation. Our mixer can operate for a maximum of 50 samples at once.

Features:

- Brushless DC motor
- Uniform acceleration
- Stepless speed regulation
- Wide speed range
- Efficient mixing

Applications:

Our Multi-tube Vortex Mixer LD-LMVM-A12 is widely used in biological research, clinical diagnostics, pharmaceuticals, and industrial applications for the efficient mixing of multiple samples.

Specifications:

Racks	Interchangeable racks
Max Load	4.5 kg
Amplitude	4 mm
Dimension	420 × 280 × 335 mm
Speed Range	500 rpm to 3000 rpm
Time Setting	1 min to 24 h
Mixing Method	Horizontal swing
Operation Mode	Touch and continuous mode
Power Consumption	100 W