



Platelet Shaker LD-LPS-A10

Labodam Equipment Ltd.

info@labodam.com | www.labodam.com

Platelet Shaker LD-LPS-A10

Overview

They are made of high-grade steel with stainless steel trays for easy cleaning and durability. It is supplied with multiple trays to accommodate 12-48 platelet bags and is low maintenance with quiet operation. They are made of high-grade steel with stainless steel trays for easy cleaning and durability. It ensures is low maintenance with quiet operation. Platelet shakers are especially designed to meet the demands of various clinical laboratories.

Features:

- Platelet bags can be loaded easily on the stainless steel basket
- Exhibits a steady and smooth 360° movement
- Attached with a tachometer for easy monitoring of speed
- Robust and practical design makes it economical choice

Applications:

Platelet shakers are widely used in blood banks and clinical research laboratories.

Specifications:

Capacity	12 bags
Power Supply	220 V, 60 Hz
Overall Dimension	370 x 255 x 390 mm
Shaking Speed Range	0 - 9 rpm