



Portable Leaf Area Meter LD-LLAM-A10

Labodam Equipment Ltd.

info@labodam.com | www.labodam.com

Portable Leaf Area Meter LD-LLAM-A10

Overview

Portable leaf area meter LD-LLAM-A10 is intended for non-destructive measuring of leaf area. Parameters like leaf area, perimeter, and aspect ratio may be measured fast and precisely. The Portable Leaf Area Meter offers a practical option for scientific research trials, plant development monitoring, and other related tasks due to its precision, efficiency, and portability.

Features:

- It may measure the area of leaves without causing any damage or altering the growth condition of plants by using a non-invasive measurement technique
- Many metrics, like leaf area, circumference, and aspect ratio, are simple to do with its user-friendly design
- Cutting-edge microprocessor technology with an LCD to visually present measurement data to increase usability
- The built-in memory can be used to store measurement data. In addition, sharing and uploading data is straightforward thanks to the USB interface

Applications:

Leaf Area Meter boasts several features to improve its usability and effectiveness for users in agriculture or plant-related research.

Specifications:

Functions	Using software and a USB connection, you may use a computer to import data, convert it to Excel format, and print it
Precision	± 2%
Test Unit	mm, cm ²
Thickness	? 8 mm
Resolution	0.01 cm ²
Gross Weight	04 kgs
Power Supply	AC 220V, 50/60 Hz
Data Capacity	? 1000 Groups
Measuring Width	? 160 mm
Measuring Length	? 1000 mm
Package Size (W x D x H)	540 x 270 x 190 mm