



Biological Microscope LD-LBM-E20

Labodam Equipment Ltd.

info@labodam.com | www.labodam.com

Biological Microscope LD-LBM-E20

Overview

Biological Microscope LD-LBM-E20 is advance research grade biological microscope with infinite corrected optical system. It offers five-field magnification power ranges from 40X to 1000X. Designed with gemel trinocular head with wide field eyepiece with diopter adjustment, Infinity Plan Semi-Apochromatic Fluorescent Objectives, DIC slide, a double layer mechanical stage, coaxial coarse and fine focusing system, N.A. 0.9 swing-out type achromatic condenser, and illumination with halogen lamp of 12 V, 100 W. Incorporated with 8.03 MP USB3.0 digital camera with advance software to capture images and videos. The camera is compatible with Windows XP/Vista/7/8/10, Mac OS X and Linux Operating System.

Features:

- Adopting color corrected infinity optical system
- Bright Field Sextuple Nosepiece, With DIC Slot for 10X, 20X, 40X objectives
- N.A.0.9 swing-out type achromatic condenser
- High-Eye point Plan Eyepiece PL10X/25 mm, Diopter Adjustable
- Infinity Plan Semi-Apochromatic Fluorescent Objectives
- Coaxial coarse and fine focusing mechanism
- Low position working stage with adjusting hand wheel

Applications:

This research grade biological microscope finds application in scientific research and testing across research units and laboratories of Medicine, Pharmacy, Life Science and other fields.

Specifications:

Eyepiece	High-Eye point Plan Eyepiece PL10X/25 mm, Diopter Adjustable
Focusing	Coaxial Focus System With Upper Limited And Tension Adjustment, Coarse Adjustment Range: 25 mm, Fine Adjustment Precision: 0.001 mm
Condenser	N.A.0.9 swing-out type achromatic condenser, with iris diaphragm, and aperture scale, provides sufficient and uniform light for full-field observation in different magnifications
Nosepiece	Bright Field Sextuple Nosepiece, With DIC Slot
Objective	Infinity Plan Semi-Apochromatic Fluorescent Objectives 4X/0.13, W.D.= 18.5 mm 10X/0.3, W.D.= 10.6 mm 20X/0.5, W.D.= 2.33 mm

	40X0.75, W.D.= 0.6 mm 100X/1.3, W.D.= 0.21 mm
Illumination	12 V, 100 W, Halogen lamp, center adjusted
Power supply	110 V to 240 V, 50 Hz/60Hz
Viewing head	Gemel trinocular head, 30° inclined, Inverted Image, Interpupillary Distance 50 to 76 mm, Three Split Ratio 0:100, 20:80, 100:0
Working Stage	Double Layer Mechanical Stage, 187×166 mm, Moving Range 80×55 mm, Precision 0.1 mm, Double Direction Transmission, Left And Right Position Handle For Option. Tension of the Torque Adjustable.
Optical system	Color corrected infinity optical system