



Distillation Tester (Front Type) LD-LDT-A14

Labodam Equipment Ltd.

info@labodam.com | www.labodam.com

Distillation Tester (Front Type) LD-LDT-A14

Overview

LDT-A14 distillation tester is a desktop structure model with adjustable heating power. It is designed according to "Test Methods for Distillation Characteristics of Petroleum Products at Barometric Pressure" as per ASTM D86.

Features:

- Desktop structure
- Flexible and convenient to operate
- Special heating furnace to ensure the safety measurement
- Adjustable heating power
- Lifting device for easy operation

Applications:

It is used to determine distillation characteristics of natural gasoline, vehicle gasoline, aviation gasoline, jet fuels, naphtha, diesel oil and similar petroleum products.

Specifications:

Dimension	620 x 510 x 620 mm
Thermometer	(-2) °C ~ 300 ° and (-2) °C ~ 400 °
Power supply	AC 220 V ± 10 %, 50 Hz
Division value	1°
Scale division	1 ml
Relative humidity	≥ 85 %
Distillation flask	125 mL
Receiving cylinder	100 ml
Ambient temperature	Room temperature to 35 °C
Heating power system	1000 W