

Dropping Point Tester (Air bath) LD-LDPT-A11

Labodam Equipment Ltd.

info@labodam.com | www.labodam.com

Dropping Point Tester (Air bath) LD-LDPT-A11

Overview

Labodam LDPT-A11 is designed and made as per the "Test Methods for Dropping Point of Lubricating Grease and Solid Hydrocarbon" • ASTM D566.

Features:

- Desktop structure
- Controlled by bidirectional silicon controlled rectifier
- Double sample working style. It can do parallel test for two samples at a time

Applications:

It is used to determine dropping point of lubricating grease and hydrocarbon over 100 °C

Specifications:

Power supply	AC 220 V ± 10 %, 50 Hz
Relative humidity	≤ 85 %
Ambient temperature	≤ 35 °C
Total power consumption	≤ 700 W
Heating power of electric furnace	250 W × 2