



**Petroleum Products Vapor Pressure Tester (REID
METHODS) LD-LVPT-C11**

Labodam Equipment Ltd.
info@labodam.com | www.labodam.com

Petroleum Products Vapor Pressure Tester (REID METHODS) LD-LVPT-C11

Overview

Labodam LVPT-C10 is designed with pressure gauge, digital display temperature controller and pump loop. It is manufactured and made as per the International Standard ASTM D323 "Test Methods for Vapor Pressure of Petroleum Products (Reid Methods)".

Features:

- Performs two experiments at same time
- Easy operation with clear display
- Stores high capacity of historical data
- Long working life
- Simple operation with elegant structure

Applications:

It is widely used in petroleum, chemical and development industry. It is suitable to determine vapor pressure of gasoline, volatile petroleum and other volatile petroleum products.

Specifications:

Test bomb	Two samples or a group of parallel samples
Power supply	AC 220 V $\pm$ 10 %
Pressure range	0 ~ 200 kPa
Bath temperature	37.8 °C
Relative humidity	$\leq$ 85 %
Ambient temperature	(-10) °C ~ 35 °C
Total power consumption	$\leq$ 1500 W
Heating power of water bath	1400 W
Temperature controlling accuracy	0.1 °C