



## **Crude Petroleum Water Content Tester LD-LWT-B10**

**Labodam Equipment Ltd.**

**[info@labodam.com](mailto:info@labodam.com) | [www.labodam.com](http://www.labodam.com)**

# Crude Petroleum Water Content Tester LD-LWT-B10

## Overview

LWT-B10 is crafted to determine water content in crude petroleum according to ASTM D4006 “Test Methods for Water Content of Crude Petroleum” international standard.

## Features:

- Simple Design
- Adjustable Electric Heating Power
- Adopts Heating Mantle

## Applications:

It is used to determine water content in crude oil, refined petroleum products, fuels, biofuels, lubricants, chemicals, and other products

## Specifications:

|                     |                        |
|---------------------|------------------------|
| Condenser           | 400 mm ± 5 mm          |
| Dimensions          | 820 × 500 × 260 mm     |
| Power Supply        | AC 220 V ± 10 %, 50 Hz |
| Water receiver      | 5 ml                   |
| scale division      | 0.05 ml                |
| Relative Humidity   | ≥ 85 %                 |
| Distillation flask  | 1000 ml                |
| Ambient Temperature | (-10) °C ~ 35 °C       |