



Pharmacy Refrigerator LD-LPRF-A16

Labodam Equipment Ltd.

info@labodam.com | www.labodam.com

Pharmacy Refrigerator LD-LPRF-A16

Overview

The Pharmacy Refrigerator LD-LPRF-A1 Series features a digital microprocessor-controlled temperature system constructed with a durable CR-Ni 304 stainless steel inner chamber, the series offers automatic defrost, efficient cooling, adjustable storage shelves, comprehensive audio-visual alarm protection, and user-friendly operation. Ideal for hospitals, pharmacies, laboratories, research institutes, and healthcare facilities.

Features:

- Uniform forced-air cooling system
- Digital temperature display
- Audible and visual alarm protection
- Lockable door with ergonomic handle
- Glass or solid insulated door options
- Adjustable storage shelves
- Energy-efficient and low-noise operation
- Optional temperature data logging and communication interface

Applications:

Widely used for storage of products like reagents, pharmaceutical compounds, laboratory grade product, medical compounds, medicinal products in pharmaceutical industries, medicals, research laboratories, etc.

Specifications:

Castors	2 Lockable + 2 Swivel
Capacity	1230 L
Door Lock	Yes
Insulation	CFC-Free Polyurethane
Refrigerant	R134a / R290
Chrome Shelf	Optional
Gross Weight	250 kg
Power Supply	220 V / 50 Hz or 110 V / 60 Hz
Control System	Microprocessor Controlled
Cooling System	Fan-Forced Air
No. of Shelves	10
PC Connectivity	USB

Thermal Printer	Optional
External Material	Electrostatic Powder-Coated Galvanized Steel
Internal Lighting	Yes
Internal Material	CR-Ni 304 Stainless Steel
Temperature Range	0°C to +15°C
Audio-Visual Alarm	Yes
Temperature Sensor	NTC
SMS & E-mail Module	Optional
Solid Insulated Door	Optional
Factory Set Temperature	+4°C
Graphical Chart Recorder	Optional
5.5-inch Touchscreen Display	Optional
Cabinet Insulation Thickness	44 mm
External Dimensions (W × D × H)	1450 × 838 × 1936 mm