



114 L Combustible Cabinet LD-LCBC-A10

Labodam Equipment Ltd.

info@labodam.com | www.labodam.com

114 L Combustible Cabinet LD-LCBC-A10

Overview

114 L Combustible Cabinet LCBC-A10 is a dual walled manual fire resistance cabinet with 38 mm of insulating air space. Doors are reinforced with a three point latch and manual lock for security. Features one adjustable shelf to accommodate storage for different containers with 2 inches vents integrated with flame arresters. Equipped with spill catcher shelves for preventing incidental drips. Inner and outer lead free powder coating provides durability and chemical resistance preventing corrosion and humidity.

Features:

- Capacity - 114 L
- Dual walled manual
- 38 mm insulating air space
- 180° fully opening door
- Three point latch and manual lock
- One adjustable shelf
- Adjustable spill retaining shelves

Applications:

Used for storage of flammable and combustible liquids away from incompatible materials (e.g., oxidizers)

Specifications:

Weight	107 kgs
Capacity	114 L
Door type	Dual door manual
Adjustable shelf	1
Dimension (H x W x D)	1120 x 1090 x 460 mm
Shelves loading capacity	50 kgs