



Corrosive Cabinet LD-LCOC-B10

Labodam Equipment Ltd.

info@labodam.com | www.labodam.com

Corrosive Cabinet LD-LCOC-B10

Overview

Corrosive Cabinet LCOC-B10 is a dual walled manual combustion cabinet with Capacity 15 L and 38 mm of insulating air space. Doors are reinforced with a three point latch and manual lock for security. Features one adjustable shelf to accommodate storage for different containers with 2 inches vents integrated with flame arresters. Equipped with spill catcher shelves for preventing incidental drips. Inner and outer lead free powder coating provides durability and chemical resistance preventing corrosion and humidity.

Features:

- Capacity - 15 L
- Dual walled manual
- 38 mm insulating air space
- 180° fully opening door
- Three point latch and manual lock
- One adjustable shelf
- Adjustable spill retaining shelves

Applications:

Used for storage of flammable and combustible liquids away from incompatible materials (e.g., oxidizers)

Specifications:

Weight	50 kgs
Capacity	15 L
Door type	Single door manual
Adjustable shelf	1
Dimension (H x W x D)	560 x 430 x 430 mm
Shelves loading capacity	50 kgs