



Corrosive Cabinet LD-LCOC-B12

Labodam Equipment Ltd.

info@labodam.com | www.labodam.com

Corrosive Cabinet LD-LCOC-B12

Overview

Corrosive Cabinet LCOC-B12 is a dual walled manual combustion cabinet with Capacity 83 L and 38 mm of insulating air space. Doors are reinforced with a three point latch and manual lock for security. Features three adjustable shelves to accommodate storage for different containers with 2 inches vents integrated with flame arresters. Equipped with spill catcher shelves for preventing incidental drips. Inner and outer lead free powder coating provides durability and chemical resistance preventing corrosion and humidity.

Features:

- Capacity - 83 L
- Dual walled manual
- 38 mm insulating air space
- 180° fully opening door
- Three point latch and manual lock
- Three adjustable shelves
- Adjustable spill retaining shelves

Applications:

Used for storage of flammable and combustible liquids away from incompatible materials (e.g., oxidizers)

Specifications:

Weight	101 kgs
Capacity	83 L
Door type	Single door manual
Adjustable shelf	3
Dimension (H x W x D)	1650 x 600 x 460 mm
Shelves loading capacity	50 kgs