



170 L Flammable Storage Cabinet LD-LFSC-A11

Labodam Equipment Ltd.

info@labodam.com | www.labodam.com

170 L Flammable Storage Cabinet LD-LFSC-A11

Overview

170 L Flammable Storage Cabinet LFSC-A11 is a dual walled manual fire resistance cabinet with 38 mm of insulating air space. Doors are reinforced with a three point latch and manual lock for security. Features two adjustable shelves to accommodate storage for different containers with 2 inches vents integrated with flame arresters. Equipped with spill catcher shelves for preventing incidental drips. Inner and outer lead free powder coating provides durability and chemical resistance preventing corrosion and humidity.

Features:

- Capacity - 170 L
- Dual walled manual
- 38 mm insulating air space
- 180° fully opening door
- Three point latch and manual lock
- Two adjustable shelf
- Adjustable spill retaining shelves

Applications:

Used for storage of flammable and combustible liquids away from incompatible materials (e.g., oxidizers)

Specifications:

Weight	148 kgs
Capacity	170 L
Door type	Dual door manual
Adjustable shelf	2
Dimension (H x W x D)	1650 x 1090 x 460 mm
Shelves loading capacity	100 kgs