



15 L Flammable Storage Cabinet LD-LFSC-B10

Labodam Equipment Ltd.

info@labodam.com | www.labodam.com

15 L Flammable Storage Cabinet LD-LFSC-B10

Overview

15 L Flammable Storage Cabinet LFSC-B10 is a dual walled manual fire resistance cabinet with 38 mm of insulating air space. Doors are reinforced with a three point latch and manual lock for security. Features one adjustable shelf to accommodate storage for different containers with 2 inches vents integrated with flame arresters. Equipped with spill catcher shelves for preventing incidental drips. Inner and outer lead free powder coating provides durability and chemical resistance preventing corrosion and humidity.

Features:

- Capacity - 15 L
- Dual walled manual
- 38 mm insulating air space
- 180° fully opening door
- Three point latch and manual lock
- One adjustable shelf
- Adjustable spill retaining shelves

Applications:

Used for storage of flammable and combustible liquids away from incompatible materials (e.g., oxidizers)

Specifications:

Weight	30.5 kgs
Capacity	15 L
Door type	Single door manual
Adjustable shelf	1
Dimension (H x W x D)	560 x 430 x 430 mm
Shelves loading capacity	50 kgs