



204 L Flammable Storage Cabinet LD-LFSC-B14

Labodam Equipment Ltd.

info@labodam.com | www.labodam.com

204 L Flammable Storage Cabinet LD-LFSC-B14

Overview

204 L Flammable Storage Cabinet LFSC-B14 is a dual walled manual fire resistance cabinet with 38 mm of insulating air space. Doors are reinforced with a three point latch and manual lock for security. Features three adjustable shelves to accommodate storage for different containers with 2 inches vents integrated with flame arresters. Equipped with spill catcher shelves for preventing incidental drips. Inner and outer lead free powder coating provides durability and chemical resistance preventing corrosion and humidity.

Features:

- Capacity - 204 L
- Dual walled manual
- 38 mm insulating air space
- 180° fully opening door
- Three point latch and manual lock
- Three adjustable shelf
- Adjustable spill retaining shelves

Applications:

Used for storage of flammable and combustible liquids away from incompatible materials (e.g., oxidizers)

Specifications:

Weight	150 kgs
Capacity	204 L
Door type	Single door manual
Adjustable shelf	3
Dimension (H x W x D)	1650 x 600 x 870 mm
Shelves loading capacity	100 kgs