



340 L Flammable Storage Cabinet LD-LFSC-C12

Labodam Equipment Ltd.

info@labodam.com | www.labodam.com

340 L Flammable Storage Cabinet LD-LFSC-C12

Overview

340 L Flammable Storage Cabinet LFSC-C12 is a dual walled self-closing fire resistance cabinet with 38 mm of insulating air space. Doors are reinforced with a three point latch and manual lock for security. Features two adjustable shelves to accommodate storage for different containers with 2 inches vents integrated with flame arresters. Equipped with spill catcher shelves for preventing incidental drips. Inner and outer lead free powder coating provides durability and chemical resistance preventing corrosion and humidity.

Features:

- Capacity - 340 L
- Dual walled self-closing
- 38 mm insulating air space
- 180° fully opening door
- Three point latch and manual lock
- Two adjustable shelves
- Adjustable spill retaining shelves

Applications:

Used for storage of flammable and combustible liquids away from incompatible materials (e.g., oxidizers)

Specifications:

Weight	186 kgs
Capacity	340 L
Door type	Dual door self-closing
Adjustable shelf	2
Dimension (H x W x D)	1650 x 11090 x 860 mm
Shelves loading capacity	100 kgs