



## **15 L Flammable Storage Cabinet LD-LFSC-D10**

**Labodam Equipment Ltd.**

**[info@labodam.com](mailto:info@labodam.com) | [www.labodam.com](http://www.labodam.com)**

# 15 L Flammable Storage Cabinet LD-LFSC-D10

## Overview

15 L Flammable Storage Cabinet LFSC-D10 is a dual walled self-closing fire resistance cabinet with 38 mm of insulating air space. Door is reinforced with a three point latch and manual lock for security. Features one adjustable shelf to accommodate storage for different containers with 2 inches vents integrated with flame arresters. Equipped with spill catcher shelf for preventing incidental drips. Inner and outer lead free powder coating provides durability and chemical resistance preventing corrosion and humidity.

## Features:

- Capacity - 15 L
- Dual walled self-closing
- 38 mm insulating air space
- 180° fully opening door
- Three point latch and manual lock
- One adjustable shelf
- Adjustable spill retaining shelf

## Applications:

Used for storage of flammable and combustible liquids away from incompatible materials (e.g., oxidizers)

## Specifications:

|                          |                          |
|--------------------------|--------------------------|
| Weight                   | 30.5 kgs                 |
| Capacity                 | 15 L                     |
| Door type                | Single door self-closing |
| Adjustable shelf         | 1                        |
| Dimension ( H x W x D )  | 560 x 430 x 430 mm       |
| Shelves loading capacity | 50 kgs                   |