



45 L Flammable Storage Cabinet LD-LFSC-D11

Labodam Equipment Ltd.

info@labodam.com | www.labodam.com

45 L Flammable Storage Cabinet LD-LFSC-D11

Overview

45 L Flammable Storage Cabinet LFSC-D11 is a dual walled self-closing fire resistance cabinet with 38 mm of insulating air space. Door is reinforced with a three point latch and manual lock for security. Features one adjustable shelf to accommodate storage for different containers with 2 inches vents integrated with flame arresters. Equipped with spill catcher shelf for preventing incidental drips. Inner and outer lead free powder coating provides durability and chemical resistance preventing corrosion and humidity.

Features:

- Capacity - 45 L
- Dual walled self-closing
- 38 mm insulating air space
- 180° fully opening door
- Three point latch and manual lock
- One adjustable shelf
- Adjustable spill retaining shelf

Applications:

Used for storage of flammable and combustible liquids away from incompatible materials (e.g., oxidizers)

Specifications:

Weight	59 kgs
Capacity	45 L
Door type	Single door self-closing
Adjustable shelf	1
Dimension (H x W x D)	890 x 590 x 460 mm
Shelves loading capacity	50 kgs