



**170 L Personal Protective Equipment
Cabinet LD-LPPE-A10**

Labodam Equipment Ltd.
info@labodam.com | www.labodam.com

170 L Personal Protective Equipment Cabinet LD-LPPE-A10

Overview

170 L Personal Protective Equipment Cabinet LPPE-A10 is a dual door storage cabinet with manual mode of operation. Equipped with three height adjustable shelves and reinforced doors. Cabinet features three-point locking system with an anti-bacterial powder coating providing durability and chemical resistance, preventing corrosion and humidity.

Features:

- Capacity - 170 L
- Reinforced doors
- Three point locking system
- Three adjustable shelves
- Anti-bacterial powder coating

Applications:

Used for the storage of personal protective equipment for emergencies, such as helmets, protective glasses, overalls etc.)

Specifications:

Weight	75 kgs
Capacity	170 L
Door type	Dual door manual
Adjustable shelf	3
Dimension (H x W x D)	1800 x 1200 x 450 mm