



114 L Weak Acid and Alkali Cabinet LD-LWAC-A10

Labodam Equipment Ltd.
info@labodam.com | www.labodam.com

114 L Weak Acid and Alkali Cabinet LD-LWAC-A10

Overview

114 L Weak Acid and Alkali Cabinet LWAC-A10 is a dual walled acid and alkali resistance cabinet with 50 mm of insulating air space. Doors are reinforced with a three point latch and manual lock for security. Features one adjustable shelf to accommodate storage for different containers. Inner and outer epoxy resin coating provides durability and chemical resistance preventing corrosion and humidity.

Features:

- Capacity - 114 L
- Dual walled manual
- 50 mm insulating air space
- Three point latch and manual lock
- One adjustable shelf

Applications:

Used to store all kinds of corrosive material or weak acid and alkali, for example: acetic acid,Carbonic acid, phosphoric acid, silicic acid, etc.

Specifications:

Weight	117 kgs
Capacity	114 L
Door type	Dual door manual
Adjustable shelf	1
Dimension (W x D x H)	1090 x 460 x 1120 mm