



170 L Weak Acid and Alkali Cabinet LD-LWAC-A11

Labodam Equipment Ltd.
info@labodam.com | www.labodam.com

170 L Weak Acid and Alkali Cabinet LD-LWAC-A11

Overview

170 L Weak Acid and Alkali Cabinet LWAC-A11 is a dual walled manual acid and alkali cabinet with 50 mm of insulating air space. Doors are reinforced with a three point latch and manual lock for security. Features two adjustable shelves to accommodate storage for different containers. Inner and outer epoxy resin coating provides durability and chemical resistance preventing corrosion and humidity.

Features:

- Capacity - 170 L
- Dual walled manual
- 50 mm insulating air space
- Three point latch and manual lock
- Two adjustable shelves

Applications:

Used to store all kinds of corrosive material or weak acid and alkali, for example: acetic acid,Carbonic acid, phosphoric acid, silicic acid, etc.

Specifications:

Weight	155 kgs
Capacity	170 L
Door type	Dual door manual
Adjustable shelf	2
Dimension (W x D x H)	1090 x 460 x 1650 mm