



340 L Weak Acid and Alkali Cabinet LD-LWAC-A13

Labodam Equipment Ltd.
info@labodam.com | www.labodam.com

340 L Weak Acid and Alkali Cabinet LD-LWAC-A13

Overview

340 L Weak Acid and Alkali Cabinet LWAC-A13 is a dual walled manual acid and alkali cabinet with 50 mm of insulating air space. Doors are reinforced with a three point latch and manual lock for security. Features two adjustable shelves to accommodate storage for different containers. Inner and outer epoxy resin coating provides durability and chemical resistance preventing corrosion and humidity.

Features:

- Capacity - 340 L
- Dual walled manual
- 50 mm insulating air space
- Three point latch and manual lock
- Two adjustable shelves

Applications:

Used to store all kinds of corrosive material or weak acid and alkali, for example: acetic acid,Carbonic acid, phosphoric acid, silicic acid, etc.

Specifications:

Weight	186 kgs
Capacity	340 L
Door type	Dual door manual
Adjustable shelf	2
Dimension (W x D x H)	1650 x 1090 x 860 mm