



45 L Weak Acid and Alkali Cabinet LD-LWAC-B11

Labodam Equipment Ltd.

info@labodam.com | www.labodam.com

45 L Weak Acid and Alkali Cabinet LD-LWAC-B11

Overview

45 L Weak Acid and Alkali Cabinet LWAC-B11 is a single walled manual acid and alkali cabinet with 50 mm of insulating air space. Doors are reinforced with a three point latch and manual lock for security. Features one adjustable shelf to accommodate storage for different containers. Inner and outer epoxy resin coating provides durability and chemical resistance preventing corrosion and humidity.

Features:

- Capacity - 45 L
- Single walled manual
- 50 mm insulating air space
- Three point latch and manual lock
- One adjustable shelf

Applications:

Used to store all kinds of corrosive material or weak acid and alkali, for example: acetic acid,Carbonic acid, phosphoric acid, silicic acid, etc.

Specifications:

Weight	66 kgs
Capacity	45 L
Door type	Single door manual
Adjustable shelf	1
Dimension (W x D x H)	700 x 570 x 1055 mm