



## **Paraffin Block Cabinet LD-LPBC-A10**

**Labodam Equipment Ltd.**  
**[info@labodam.com](mailto:info@labodam.com) | [www.labodam.com](http://www.labodam.com)**

# Paraffin Block Cabinet LD-LPBC-A10

## Overview

Paraffin Block Cabinet is metallic standalone unit, comes with 18 drawers. Paraffin cabinet is used to store paraffin blocks, Drawers is consisting of ABS slide plate and cabinet material is cold rolled steel plate and surface is treated with Degree oil, Table regulator, Zinc Parkerizing, Passivation, Powder spraying which is resistant to most chemical aggression.

## Features:

- 18 Drawers with loading capacity of 12000 paraffin block
- Slot is made by special mold
- Corrosion resistant material

## Applications:

Paraffin block cabinet is used for efficient and organized storage of paraffin blocks across histology, pathology, forensics, research centers etc.

## Specifications:

|                           |                                                                                                                         |
|---------------------------|-------------------------------------------------------------------------------------------------------------------------|
| Capacity                  | 18 Drawers (3 layers x 6 grid, plus a base)<br>Loading: 12000 pcs paraffin block                                        |
| Mark groove               | Integration forming                                                                                                     |
| Gross weight              | 29kg, 3pcs; 33kg, 1pc                                                                                                   |
| Slot material             | Made by special mold                                                                                                    |
| Drawer material           | Specialized slide drawer, ABS slide plate                                                                               |
| Cabinet material          | Cold-rolled steel plate. Surface Treatment: Degree oil, Table regulator, Zinc Parkerizing, Passivation, Powder spraying |
| Package Size (W x D x H)  | 530 x 570 x 440 mm, 3pcs; 530 x 570 x 565 mm, 1pc                                                                       |
| External Size (W x D x H) | 450 x 478 x 1295 mm                                                                                                     |