



Slide Storage Cabinet LD-LSLC-A10

Labodam Equipment Ltd.

info@labodam.com | www.labodam.com

Slide Storage Cabinet LD-LSLC-A10

Overview

Slide Storage Cabinet LD-LSLC-A10 is metallic standalone unit, comes with 18 drawers. Slide cabinet is used to store glass slides, Drawers is consisting of ABS slide plate and cabinet material is cold rolled steel plate and surface is treated with Degree oil, Table regulator, Zinc Parkerizing, Passivation, Powder spraying which is resistant to most chemical aggression.

Features:

- 4563 or 9072 glass slides can be stored
- 18 Drawers Cabinet
- Slot is made by special mold
- Corrosion resistant material
- Large storage capacity

Applications:

Slide storage cabinet is used for efficient and organized storage of glass slides across histology, pathology, forensics, research centers etc.

Specifications:

Capacity	18 Drawers (3 layers x 6 grid, plus a base) Loading: 4536 or 9072 pcs glass slides
Mark groove	Integration forming
Gross weight	27 kgs, 3pcs; 31 kgs, 1pc
Slot material	Made by special mold
Drawer material	Specialized slide drawer, ABS slide plate
Cabinet material	Cold-rolled steel plate
Package Size (W x D x H)	530 x 570 x 440 mm, 3pcs; 530 x 570 x 565mm, 1pc
External Size (W x D x H)	450 x 478 x 1295 mm