



COBB sample cutter LD-TP-C30

Labodam Equipment Ltd.

info@labodam.com | www.labodam.com

COBB sample cutter LD-TP-C30

Overview

COBB sample cutter TP-C30 with user friendly operations, benchtop design, is an important instrument for paper & cardboard cutting. It can cut the sample quickly & precisely. It is the ideal auxiliary instrument for the packaging, paper industry.

Features:

- Circular sample cutter
- Maximum sample diameter 125 mm
- Sample thickness range 0.1 to 1.0 mm
- User friendly operations

Applications:

COBB sample cutter is used in paper & packaging industry for quick and precise cutting of paper, board.

Specifications:

Dimension	240×288×435 mm
Net Weight	25 kg
Sample error	± 0.2 mm
Sample diameter	125 mm
Sample Thickness	0.1 to 1.0 mm