



Coefficient of friction tester LD-TP-D10

Labodam Equipment Ltd.
info@labodam.com | www.labodam.com

Coefficient of friction tester LD-TP-D10

Overview

Coefficient of friction tester TP-D10 is used in testing the static & kinetic friction coefficient of plastic films and sheet. Tester is designed in accordance with ISO 8295 standard. It is supported by professional software to print single or group results, reports and curves automatically.

Features:

- Microprocessor controlled computer system
- Equipped with build-in mini printer
- Quick start options for starting new test
- Simultaneously measures kinetic & static coefficients of friction
- LCD screen display
- Easy to use

Applications:

Coefficient of friction tester is used to test kinetic and static coefficient of friction samples, materials used in paper & cardboard industry.

Specifications:

Slider size	63×63 mm
Slider stroke	Max 100 mm (adjustable)
Slider weight	200 g
Speed of slider	100 ± 1 0mm/min
Testing accuracy	± 2 %
Measurement range	0.01 to 10 N
Dimension of working plate	170×225 mm