



## **Standard strip cutter LD-TP-G10**

**Labodam Equipment Ltd.**

**[info@labodam.com](mailto:info@labodam.com) | [www.labodam.com](http://www.labodam.com)**

# Standard strip cutter LD-TP-G10

## Overview

Standard strip cutter TP-G10 with sleek and advanced design is used to cut the sheet of cardboard or paper into fine strips. It's a common standard cutter which can cut the sample with 15 mm width. With user-friendly operation, it finds applications in paper & printing industry.

## Features:

- Advanced and sleek design
- Sample width 15 mm
- Sample maximum length 300 mm
- Easy to use

## Applications:

Strip cutter is used in paper, packaging, inspection department, for paper crafts making industry.

## Specifications:

|                       |               |
|-----------------------|---------------|
| Weight                | 4 kg          |
| Dimension             | 240×520×70 mm |
| Parallelism           | ± 0.1 mm      |
| Sample width          | 15 mm         |
| Sample maximum length | 300 mm        |