



Notch Cutter LD-TNC-A10

Labodam Equipment Ltd.

info@labodam.com | www.labodam.com

Notch Cutter LD-TNC-A10

Overview

Notch Cutter TNC-A10 is used in sample preparation of sample for Izod & charpy impact testing of the plastic & other non- metal material. It is a superior & precise sample preparation machine with simple structural body & convenient operations. It is in accordance with international standards such as ISO 179, ASTM D256.

Features:

- Cutting tool distance: 16 mm
- Can make notch sample for Izod as well as charpy test
- Samples setting with different tool size
- User friendly design

Applications:

It is use for scientific and research work across colleges and institutes, across factories for production etc. to prepare notch making samples.

Specifications:

Gross weight	33 kg
Dimension (L*W*H)	500×200×300 mm
Cutting tool distance	16 mm
Moving forward distance	More than 20 mm
Cutting tool parameter Model A	$45^{\circ} \pm 30'$ & $r=0.25 \pm 0.025$
Cutting tool parameter Model B	$45^{\circ} \pm 30'$ & $r=1.0 \pm 0.025$
Cutting tool thickness Model C	$H=2 \pm 0.1$ & $r \geq 0.1$