



Multifunctional Ultrasonic Cleaner LD-LMFU-A14

Labodam Equipment Ltd.

info@labodam.com | www.labodam.com

Multifunctional Ultrasonic Cleaner LD-LMFU-A14

Overview

Multifunctional Ultrasonic Cleaner LMFU-A14 is a single 30 L multifunctional tank with heat and ultrasonic action providing deep cleansing, removing contamination even from blind holes and internal surfaces. Its multifunctionality includes power adjusting mode, pulse or degas control and variable ultrasonic intensity control. Equipped with degassing and degreasing function for removal of dissolved gasses and cleansing of the tank.

Features:

- 30 L capacity tank
- Integrated digital circuit for timer and heater
- Temperature range : 0 °C ~ 100 °C
- Power adjustable mode
- Cleaning duration : 1 ~ 99 mins (can be kept "ON" always)
- 4 working modes available (Power , Pulse , Turbo , Degas)
- One key access for auto stand-by , sleep and wake up mode
- On power cut, auto memory save option
- Efficient cleaning by distilled water
- Equipped with degassing and degreasing function
- Precise digital processed temperature controls
- Optimum durability under prolonged and heavy usage

Applications:

Ideal for cleaning irregular and complex instruments, lenses and other optical parts, dental and surgical instruments, electronics. It is also used in extraction and to degas certain analyzed stuff in laboratories

Specifications:

Timer	1 ~ 99 mins (optional : always ON)
Voltage	AC 220 V / 50 Hz
Capacity	30 L
Drainage	Yes
Heater power	500 W
Temperature set	0 °C ~ 100 °C
Ultrasonic power	0 ~ 600 W
Construction material	Stainless Steel
Tank size (L x W x H)	500 x 300 x 200 mm
Overall size (L x W x H)	550 x 330 x 407 mm