



Power Adjustable Ultrasonic Cleaner LD- LPAU-A14

Labodam Equipment Ltd.
info@labodam.com | www.labodam.com

Power Adjustable Ultrasonic Cleaner LD-LPAU-A14

Overview

Power Adjustable Ultrasonic Cleaner LPAU-A14 is the most reliable instrument with 30 L of tank capacity and cleaning duration of 1 ~ 99 mins for efficient intensive cleaning of instruments and products making them contamination, rust and germ free. Its key feature is the variance of using it at 28 KHz for gross rinse to 40 KHz for intensive rinse. The power can be adjusted which provides the user with a precise control from zero watts to maximum power.

Features:

- 30 L capacity tank
- Temperature range : 0 °C ~ 80 °C
- Cleaning duration : 1 ~ 99 mins
- 28 KHz for gross rinse to 40 KHz for intensive rinse
- One key access for auto stand-by , sleep and wake up mode
- Integrated digital circuit for timer and heater
- Precise digital processed temperature controls
- Efficient cleaning by distilled water
- Optimum durability under prolonged and heavy usage
- Equipped with degassing and degreasing function

Applications:

Ideal for cleaning irregular and complex instruments, lenses and other optical parts, dental and surgical instruments, electronics. It is also used in extraction and to degas certain analyzed stuff in laboratories

Specifications:

Timer	1 ~ 99 mins
Voltage	AC 220 V / 50 Hz
Capacity	30 L
Drainage	Yes
Heater power	500 W
Temperature set	0 °C ~ 80 °C
Ultrasonic power	0 ~ 600 W
Construction material	Stainless Steel
Tank size (L x W x H)	500 x 300 x 200 mm
Overall size (L x W x H)	550 x 330 x 387 mm