



Mini Vacuum Pump LD-LMVP-A10

Labodam Equipment Ltd.

info@labodam.com | www.labodam.com

Mini Vacuum Pump LD-LMVP-A10

Overview

Mini vacuum pump LMVP-A10 uses big preface torque configuration and is particularly designed for low temperature working environment. It portable and requires less storage place. This mini vacuum pump has space-saving design and is ideal for systems integration.

Features:

- Compact, reliable and powerful
- Especially designed for low temperature working environment
- High pumping speeds with small dimensions
- Eco-friendly exhaust gas free of oil mists
- Quiet, low-vibration operation
- Easy to maintain

Applications:

It can be used for vacuum packaging, vacuum freezing, gas analysis and other vacuum evacuations.

Specifications:

Power	90 / 120 W
Weight	6.5kg
Dimension	243 x 114 x 207 mm
Oil Capacity	220 ml
Power Supply	110 V, 60 Hz / 220 V, 50 Hz
Pumping Speed	28.31 to 22.65 L / min
Rotating Speed	1720 / 1440 rpm
Ultimate Pressure	5 Pa