

Benchtop Conductivity Meter LD-LBCM-A15

Labodam Equipment Ltd.

info@labodam.com | www.labodam.com

Benchtop Conductivity Meter LD-LBCM-A15

Overview

Benchtop Conductivity Meter LD-LBCM-A15 with 1 to 3 points push button conductivity calibration, automatically recognizes calibration standards. With setup menu, it does setting of 6 parameters calibration points, temperature units and others. 7 inch TFT display shows settings and measurement values. Temperature compensation by manual & automatic method improves the accuracy of conductivity meter.

Features:

- **Common Features**
- Stability indicator displays the status of current measurement
- Measures Conductivity, TDS, Salinity, Resistivity and Temperature of the liquids
- Automatic temperature compensation for accurate values over the entire range
- Selectable temperature units (°C or °F) for different application requirements
- Reset function automatically resumes all the setting back to factory default options
- Manual and automatic mode of temperature calibration
- Auto hold feature to lock, record and read stable measurement endpoint
- Password protection to avoid unauthorized access to settings
- Expanded memory stores and recalls up to 1000 readings
- **Conductivity / TDS / Salinity / Resistivity Mode**
- 1 to 3 Conductivity push button calibration points
- 1 Temperature calibration point
- Can select cell constant (0.1/1/10) and normalization temperature (20°C or 25°C)
- TDS conversion factor, seawater, practical salinity measurements mode to meet different water quality measurement requirements
- Linear or pure water compensation improves the measurement accuracy

Applications:

Used for quality checking of water on different parameters across research, QC labs, water management plants and others production units.

Specifications:

Memory	Stores Up to 1000 data sets						
Output	USB interface						
Weight	1.7 kg						
Display	LCD						
TDS Mode	<table><tr><td>TDS measurement Range</td><td>0 to 10.00, 100.0, 1000 mg/L, 10.00, 200.0 g/L</td></tr><tr><td>TDS Accuracy</td><td>± 1 % F.S.</td></tr><tr><td>TDS Factor</td><td>0.01 to 1.00</td></tr></table>	TDS measurement Range	0 to 10.00, 100.0, 1000 mg/L, 10.00, 200.0 g/L	TDS Accuracy	± 1 % F.S.	TDS Factor	0.01 to 1.00
TDS measurement Range	0 to 10.00, 100.0, 1000 mg/L, 10.00, 200.0 g/L						
TDS Accuracy	± 1 % F.S.						
TDS Factor	0.01 to 1.00						

	(default 0.5)</td></tr><tr><td>TDS Resolution</td><td>0.01 / 0.1 / 1 or automatic</td></tr></table>												
Connector	6-pin mini plug												
Power supply	DC 12V, Power Adapter												
Salinity Mode	<table> <tr> <td>Salinity Range</td> <td>0.00 to 80.00 ppt / 0.00 to 42.00 psu / 0.00 to 8.00%</td> </tr> <tr> <td>Salinity Resolution</td> <td>0.01</td> </tr> <tr> <td>Salinity Accuracy</td> <td>± 1 % F.S.</td> </tr> </table>	Salinity Range	0.00 to 80.00 ppt / 0.00 to 42.00 psu / 0.00 to 8.00%	Salinity Resolution	0.01	Salinity Accuracy	± 1 % F.S.						
Salinity Range	0.00 to 80.00 ppt / 0.00 to 42.00 psu / 0.00 to 8.00%												
Salinity Resolution	0.01												
Salinity Accuracy	± 1 % F.S.												
Reset Function	Yes												
Resistivity Mode	<table> <tr> <td>Resistivity Range</td> <td>0.00 to 30.00 M?</td> </tr> <tr> <td>Resistivity Resolution</td> <td>0.01 / 0.1 / 1 or automatic</td> </tr> <tr> <td>Resistivity Accuracy</td> <td>± 1 % F.S.</td> </tr> </table>	Resistivity Range	0.00 to 30.00 M?	Resistivity Resolution	0.01 / 0.1 / 1 or automatic	Resistivity Accuracy	± 1 % F.S.						
Resistivity Range	0.00 to 30.00 M?												
Resistivity Resolution	0.01 / 0.1 / 1 or automatic												
Resistivity Accuracy	± 1 % F.S.												
Temperature Mode	<table> <tr> <td>Temperature range</td> <td>0 to 105°C, 32 to 221°F</td> </tr> <tr> <td>Temperature accuracy</td> <td>± 0.5 °C</td> </tr> <tr> <td>Temperature resolution</td> <td>± 1 °C</td> </tr> <tr> <td>Temperature calibration points</td> <td>1 point</td> </tr> <tr> <td>Temperature coefficient</td> <td>0.0 to 10.0 % per °C</td> </tr> <tr> <td>Temperature Compensation Range</td> <td>0 to 100°C / 32 to 212°F</td> </tr> </table>	Temperature range	0 to 105°C, 32 to 221°F	Temperature accuracy	± 0.5 °C	Temperature resolution	± 1 °C	Temperature calibration points	1 point	Temperature coefficient	0.0 to 10.0 % per °C	Temperature Compensation Range	0 to 100°C / 32 to 212°F
Temperature range	0 to 105°C, 32 to 221°F												
Temperature accuracy	± 0.5 °C												
Temperature resolution	± 1 °C												
Temperature calibration points	1 point												
Temperature coefficient	0.0 to 10.0 % per °C												
Temperature Compensation Range	0 to 100°C / 32 to 212°F												
Conductivity Mode	<table> <tr> <td>Conductivity Range</td> <td>0.01 to 20.00, 200.0, 2000 µS/cm, 20.00, 200.0 mS/cm</td> </tr> <tr> <td>Conductivity Accuracy</td> <td>±0.5 % F.S.</td> </tr> <tr> <td>Conductivity Resolution</td> <td>0.001/0.01/0.1/1, automatic</td> </tr> <tr> <td>Conductivity Calibration Points</td> <td>1 to 3 Points</td> </tr> <tr> <td>Conductivity Solutions</td> <td>10, 84, 1413 µS/cm, 12.88,</td> </tr> </table>	Conductivity Range	0.01 to 20.00, 200.0, 2000 µS/cm, 20.00, 200.0 mS/cm	Conductivity Accuracy	±0.5 % F.S.	Conductivity Resolution	0.001/0.01/0.1/1, automatic	Conductivity Calibration Points	1 to 3 Points	Conductivity Solutions	10, 84, 1413 µS/cm, 12.88,		
Conductivity Range	0.01 to 20.00, 200.0, 2000 µS/cm, 20.00, 200.0 mS/cm												
Conductivity Accuracy	±0.5 % F.S.												
Conductivity Resolution	0.001/0.01/0.1/1, automatic												
Conductivity Calibration Points	1 to 3 Points												
Conductivity Solutions	10, 84, 1413 µS/cm, 12.88,												

	111.8 mS/ cm</td></tr><tr><td>Conductivity Temperature coefficient</td><td>0.0 to 10.0 %/°C</td></tr><tr><td>Compensation coefficient</td><td>Linear / Non-Liner / Pure Water</td></tr><tr><td>Cell Constant</td><td>2-Pole Electrode K=0.1, 1, 10 or 4-Pole Electrode</td></tr><tr><td>Reference temperature</td><td>20 or 25°C</td></tr></table>
Data Hold Function	Manual or Automatic Endpoint Detection
Stability Criteria	Standard or High
Calibration Due Alarm	1 to 31 Days
Dimension (L × W × H)	240 × 220 × 80 mm