



Adherent Cell Culture Flask LD-ACF104L

Labodam Equipment Ltd.

info@labodam.com | www.labodam.com

Adherent Cell Culture Flask LD-ACF104L

Features:

- Material: High Transparency Polystyrene
- Tissue culture treated
- Sterilization Level is 10^{-6}
- Available in both vent cap and plug seal cap
- Gamma ray sterilized
- DNase, RNase and pyrogenic free

Specifications:

Cap	Plug
Capacity	600 mL
Culture area	175 cm^2
Quantity per case	40