



Adherent Culture Plate LD-ACP100L

Labodam Equipment Ltd.

info@labodam.com | www.labodam.com

Adherent Culture Plate LD-ACP100L

Features:

- Material: Polystyrene
- Sterilization level is 10^{-6}
- Tissue Culture Treated
- Number codes for easy identification
- Flat surface bottom
- Maximized growth area
- Raised rims on wells to avoid cross contamination
- DNase, RNase and pyrogenic free

Specifications:

Packaging	50 pcs/blister
No of Well	6
Culture area	9.61 cm^2
Packaging Type	Blister
Quantity per case	200