



Adherent Culture Plate LD-ACP107L

Labodam Equipment Ltd.

info@labodam.com | www.labodam.com

Adherent Culture Plate LD-ACP107L

Features:

- Material: Polystyrene
- Sterilization level is 10^{-6}
- Tissue Culture Treated
- Number codes for easy identification
- Flat surface bottom
- Maximized growth area
- Raised rims on wells to avoid cross contamination
- DNase, RNase and pyrogenic free

Specifications:

Packaging	50 pcs/bag
No of Well	96
Culture area	0.35 cm^2
Packaging Type	Bag
Quantity per case	200