



## **Polyoxymethylene Embedding Cassette LD-PEC103L**

**Labodam Equipment Ltd.**

**[info@labodam.com](mailto:info@labodam.com) | [www.labodam.com](http://www.labodam.com)**

# Polyoxymethylene Embedding Cassette LD-PEC103L

## Features:

- Material: Polyoxymethylene
- Cassette writing face is at a 30° and 45° angle
- Cover and base can be disassembled and reassembled with ease
- Resistant to solvents

## Specifications:

|                   |                           |
|-------------------|---------------------------|
| Type              | Strip Holes               |
| Color             | Green                     |
| Lid Type          | Detachable Lid            |
| Packaging         | 250 pcs/Bag, 20 Bags/Case |
| Quantity Per Case | 5000                      |