



Polyoxymethylene Embedding Cassette LD-PEC105L

Labodam Equipment Ltd.

info@labodam.com | www.labodam.com

Polyoxymethylene Embedding Cassette LD-PEC105L

Features:

- Material: Polyoxymethylene
- Cassette writing face is at a 30° and 45° angle
- Cover and base can be disassembled and reassembled with ease
- Resistant to solvents

Specifications:

Type	Fine Square Holes
Color	White
Lid Type	Detachable Lid
Packaging	250 pcs/Bag, 20 Bags/Case
Quantity Per Case	5000